




## Alak A.V. Fistula Needle



### Features

- Ultra wall and siliconized fistula needle with or without back eye.
- The stainless steel needle with ultra thin wall
- Smooth kink resistant tube.
- For easy gripping & fixation comes with flexible wings
- Color wing for easy identification of needle size.
- Sterile / Single Use/ Individual Packed

Colour Coded	Gauge	Needle	
		Ext.Dia.mm	Length mm

Blue	G-15	1.84	38
Green	G-16	1.65	38
Orange	G-17	1.47	38

**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com