




**ALAK**  
HEALTHCARE PVT. LTD.

## Alak Infusion Set (I.V.SET)





## Features

- Clear, Transparent & Flexible drip Chamber rapid adjustment of fluid
- Extra Smooth kink resistance PVC tube 150 cm
- Disc Filter in the drip chamber.
- Flow rate =20 drops/ml +/- 0.1ml
- Very Efficient roller controller for accurate and unrestricted flow
- Available with a Y-injection site (with or without Latex)
- As per International standard 6% taper Luer.
- Option of Luer locking , Luer slip , with & without Airvent , latex bulb , latex tube , Needle any size G : 21
- Sterile / Disposable / Individual Packed , Non-Pyrogenic & Non Toxic
- Accurate flow controller.

## Packing

Polyethylene pack, Paper pack, Ribbon pack, Box pack available


**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



## Alak Infusion Set (I.V.SET)



### Features

- Clear, Transparent & Flexible drip Chamber rapid adjustment of fluid
- Extra Smooth kink resistance PVC tube 150 cm
- Disc Filter in the drip chamber.
- Flow rate =60 drops/ml +/- 0.1ml
- Very Efficient roller controller for accurate and unrestricted flow
- Available with a Y-injection site (with or without Latex)
- As per International standard 6% taper Luer.
- With Luer locking , with & without Airvent , latex bulb , latex tube , Needle any size G :23
- Sterile / Disposable / Individual Packed , Non-Pyrogenic & Non Toxic
- Accurate flow controller

### Packing

Polyethylene pack, Paper pack, Ribbon pack, Box pack available


**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



REF NO.	DESCRIPTION OF GOODS	PACKING SIZE INNER BAG
AL IV-001	Infusuin Set (IV Sets) Without Airvent with latex tube & needle (Economy Pack)	25 PCS
AL IV-002	Infusuin Set (IV Sets) With Airvent (vented spike) with latex bulb & needle	50 PCS
AL IV-003	Infusuin Set (IV Sets) Without Airvent with latex bulb & needle	50 PCS
AL IV-004	Infusuin Set (IV Sets) With or Without Airvent with Y injection site (latex) & Luer Lock	50 PCS
AL IV-005	Infusuin Set (IV Sets) With or Without Airvent with Needle Free Y site & Luer Lock	50 PCS
AL IV-006	Infusuin Set (IV Sets) With Flow Regulator With Y site (latex) & Luer Lock	50 PCS
AL IV-007	Infusuin Set (IV Sets) With & Without Airvent with 3 way stop cock & Luer Lock	50 PCS



**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com