




## Alak Nelaton Catheter



### Details

- Nelaton catheters are used for short term bladder catheterisation.
- Manufactured from non-toxic, non-irritant medical grade soft P.V.C.
- Compatible with catheter lubricants.
- Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage
- Available with the option of Radio opaque line throughout the length for X-ray visualization
- Frozen surface tubing for super smooth intubation.
- Proximal end is fitted with universal funnel shaped connector for extension.
- Funnel shape connectors are color coded for instant identification of sizes.
- Available in Male and Female versions.
- Sterile product for single use
- 38cm Nelaton catheter / 23cm female catheter

Color Coded									
	Blue	Black	White	Green	Orange	Red	Yellow	Violet	Royal Blue
Size	8	10	12	14	16	18	20	22	24

**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.


**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



**ALAK**  
HEALTHCARE

FG



**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com