



Alak Scalp Vein Set



Features

- Rapid venous access, for patient comfort during infusion.
- Short beveled siliconized needle facilitates traumatic cannulation
- Thin wall high quality needle provides better flow rate per gauge,
- Butterfly wings facilitate for easy fixation with the Skin
- Soft non –toxic non irritant tube which does not kink or coil set
- Proximal end is provided with flexible female luer fitting
- Supper smooth kink resistant tube.
- Color coded butterfly wing by size for easy identification.

Size	G-18	G-19	G-20	G-21	G-22	G-23	G-24	G-25
Colour Code	Pink	Brown	Yellow	Green	Black	Sea Green	Red	Blue

Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com