




Alak Twin Bore Nasal Oxygen Set



Details

- Manufactured from Soft, non toxic PVC Thus non irritating even in long term use.
- Soft twin prong nasal tips to ensure equal volume of oxygen to both air passages.
- Smoothly finished & adjustable nasal tips for maximum patient's comfort.
- Star lumen tube ensures the supply of oxygen even if tube is accidently kinked.
- Soft Funnel shaped connector facilitates easy connection to oxygen sources.
- Sterile/ Disposable/ Individual/ Packed.