



Alak Male External Catheter



Widely used by urologists and surgeons

- Made of soft, comfortable rubber.
- For urine incontinence for day and night use in male patient.
- Provided with self adhesive coated strip.
- Designed to connect to Urine Collection Bags without internal catheterization.
- Sterile / Disposable / Individual Packed.

Sizes: Extra Small, Small, Medium, Large and Extra Large

Small	Medium	Large	Ex. Large
20 mm	25 mm	30 mm	35 mm



Alak Female Catheter



- Manufactured from soft none toxic, non irritant kink resistant medical grade PVC tube.
- Fitted with Universal Funnel shaped connector for safe connection to Urine bag.
- Color coded for easy identification of size.
- Option of Radio opaque line throughout the length for X-ray visualization.
- Closed tip with two lateral eyes for efficient drainage.
- Sterile / Single Use / individual packed.









Color Coded								
	Light Green	Blue	Black	White	Green	Orange	Red	Yellow
Size FG	6	8	10	12	14	16	18	20



Alak Rectal Catheter



- Made of none toxic, non irritant kink resistant medical grade PVC tube.
- Funnel shaped connector for easy connection to the suction equipment.
- Option of Radio opaque line throughout the length for X-ray visualization.
- Closed tip with two lateral eyes for efficient drainage.
- Color coded for easy identification of size.
- Sterile / Disposable / Individual packed.

								
	Red	Yellow	Violet	Light Blue	White	Green	Light Grey	Brown
	18	20	22	24	26	28	30	32