



Alak Chest Drainage Bottle



Features

- Designed to forward the collected Drainage into a collection chamber by the liquid static suction effects.
- Easy to read graduation for drain volume
- Facilitates safe Retrograde Technique
- Easy to observe and measure the collective drainage content during surgery.
- The hanging at a higher position as compares to the common collection chamber allowing the collection of drainage without any inconvenience.
- Suction port is provided to fill water.
- Kink resistant large bore tubing facilities.
- Unrestricted flow.

Also available The Kit for Safe mobilisation of patients when under water seals cannot be established.

- Ambulatory chest drainage system with integral tubing.
- Syringe for primary flutter valve.
- Inner sheet sterile field.