




ALAK
HEALTHCARE

Clips & Appliers

Titanium Clips



The Clip Applier is a reusable instrument designed to deliver and close metallic legating clips. This instrument is combined with Hemosec Clips. Is made from either 316 L stainless steel or Titanium. The legating clips are positioned around a tubular structure and closed by applying pressure to the handles of the applier.


Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



Open Applier



The applier has ergonomically designed handle that provides comfort and control. It has diamond shaped grooving for a traumatic secure ligation
Triangular shaped wire for easy squeeze application.



ALAK
HEALTHCARE

Laparoscopic Applier



Laparoscopic Applier which is a 35cm long disposable cannulated applicator in a pack of 6. Device obturator allows for a precise FloSeal application to the base of the lesion. Handle ergonomics facilitate device manipulation. Cannula rigidity aids in device navigation and precise material delivery.

Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com