




Alak Endotracheal Tube (Cuffed / Plain)



Endotracheal Tube (Cuffed)

Details

- Manufactured from siliconised non-toxic, non-irritant P.V.C. kink resistant thermo sensitive tube.
- Kink resistant inflation tube ensures patient safety during cuff inflation and deflation.
- Soft low pressure cuff provides an air tight fit between tube and trachea.
- X-ray opaque line.
- Graduated marking to display distance from the tip.
- Standard 15mm connector provided at proximal end.

Plain

Size	2	2.5	3	3.5	4	4.5	5	5.5
------	---	-----	---	-----	---	-----	---	-----

Cuffed

Size	4	4.5	5	5.5	6	6.5
7	7.5	8	8.5	9	9.5	10