



Alak Infant Feeding Tube



Features

- For neonatal and pediatric nutritional feeding
- Distal end is coned with two lateral eyes.
- Proximal end is provided with female luer mount.
- Low friction tubing and super smooth tip in straight pack ensure trauma free intubation.
- Flexible connector with cap to prevent reflux and contamination
- For Identification of the respective size the connectors of the knotted are of different colors.
- Radio- opaque line is provided throughout its length for X-ray visualization.
- Manufactured from non-toxic P.V.C.
- Non irritant to delicate mucosa.
- Tube Length 52 cm
- Size : 5/6/7/8/9/10 FG
- Sterile / Disposable / individual packed.

Color Coded						
	Grey	Green	White	Blue	Orange	Black
Size FG	5	6	7	8	9	10