



**ALAK**  
HEALTHCARE

**ALAK IV Cannula / Catheter**  
**With Wings & Injection Port**



**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



## Features

- Needle: Siliconised stainless steel back cut needle to perform easy and triple point for painless insertion.
- Catheter: Made of Biocompatible FEP material has smooth outer surface and tapered tip for very smooth insertion with minimum trauma.
- The Flash Back chamber allows quick visualization of venous return.
- Flexible Wings: The flexible wings for easy gripping and safe clamping.
- Needle cover: Soft and transparent needle cover to prevent accidental damage to needle & catheter
- Sterile , Single Use & Non Pyrogenic
- Options Available:
  - FEP/PTFE/PUR Catheter
  - Hydrophobic filter
  - Suturable Wings
  - Color coded body / port cap as per ISO Standard.
- Hard Blister pack using Medical Grade paper / Typek with standard Plastic Film. Soft Blister pack using with PVC Free Plastic Film.

SPECIFICATIONS				
Gauge	Colour Code	Ext. Dia. mm	Length mm	Flow Rate ml/min
14G	Orange	2.1	45	270
16G	Grey	1.8	45	180
18G	Green	1.3	45	90
20G	Pink	1.1	32	65
22G	Blue	0.9	25	36
24G	Yellow	0.7	19	23
26G	Violet	0.6	19	15