



Alak Levin's Tube



Features

- Suitable for gastro intestinal feeding and aspiration.
 - Distal end is coned with four lateral eyes.
 - Radio opaque line throughout its length helps checking the position radio logically.
 - Low friction medical grade P.V.C. reduces the incidence of nasal discharge.
 - Proximal end is provided with universal funnel.
 - Marked at 45, 55, 65 & 75 cm from the distal end for accurate placement
- Sterile / Single Use / Pyrogen Free / Individual packed.

Sizes : FG 8/10/12/14/16/18/20 with Color Code

| | | | | | | | |
|-------------|------|-------|-------|-------|--------|-----|--------|
| Color Coded | | | | | | | |
| | Blue | Black | White | Green | Orange | Red | Yellow |
| Size FG | 8 | 10 | 12 | 14 | 16 | 18 | 20 |