



ALAK Low & High Pressure Extension Tube



Designed for transfusion and infusion to extend the IV fluid pathway & for Blood / Blood derivatives

Features

- Sterile, Disposable and Non- Pyrogenic
- Two pressure variants 60 PSI & 100 PSI.
- 6% luer taper.
- With or Without "Y" Connection.
- Male luer lock at one end and female luer lock at other end.
- Soft and kink resistant PVC tubing
- 3.0*4.1mm tube with constant dimensions.
- Tube Length: 25, 50, 100, 150, & 200 cm.

EXTENSION TUBE with T Connector with Tube Length: 15 cm