




## Alak Measured Volume Set

### Burette Set (110ml / 150 ml)



### Features

- Flexible, Transparent & soft cylindrical measured volume chamber of 110 ml / 150ml
- Beveled strong spike for smooth penetration.
- Floating auto shut off valve acts as floating indicator to prevent air trapping in the fluid line.
- Available with or without "Y" connection , Bulb Latex , Tube Latex , Latex Free , Luer Tip / Luer Lock
- Soft and kink Resistant tube.
- Side Clamp / Pinch Clamp for extra medication and continuous change over.
- With or Without Needle , Size G: 23.
- Eto Sterile / Disposable / Individual Packed.
- Drop size of 60 drops/ml
- This device is useful for giving saline to the infants or heart patient



MEASURED VOLUME ADMINISTRATION SETS ( BURETTE SETS)		
Ref No.	Description Of Goods	Packing Size Inner Bag
AL BU-101	Measured Volume Administration Set ( Burette ) Size 100 ml/ 110ml/ 150 ml vented / nonvented spike, latex or latex free, with luer lock & needle	10 PCS
AL BU-102	Measured Volume Administration Set ( Burette ) Size 100 ml/ 110ml/ 150 ml vented / nonvented spike, With Y site luer lock & needle	10 PCS
AL BU-103	Measured Volume Administration Set ( Burette ) Size 100 ml/ 110ml/ 150 ml vented / nonvented spike, With Y site & luer lock & with flow regulator	10 PCS

