



ALAK Nebulizer / Breath Easy



Details

- Compressor based Nebulizer machine to administer medication in the form of a mist directly to the patient's airways.
 - Allow the medication to go directly to a person's airway and lungs
 - Inhaling medication easier for people who are having serious difficulty in breathing. or who have trouble using an inhaler which requires careful timing and control of one's breathing.
1. Compressor based
 2. Compact
 3. Useful for all ages
 4. Can be operated continuously for 20 min- ON & after 40min-OFF

Specification:

Maximum Compressor Pressure	: 2.5-3.0 bar
Nebulizer Rate	: 0.3 ml per min
Nebulizer Cup Capacity	: 6 ml
Particle Size	: 0.5 to 5 (miu) m
Noise level	: 55 dba
Weight (compressor only)	: 1.83 kg (4.0 pound)
Electrical Requirement	: 120v/60 h2-220v/ 50h2

Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



Nebulizer Assembly & Accessories:

- Nebulizer cup
- Tubing
- An adult mask Mouthpiece
- Nebulizer insert
- Nebulizer cap
- Extra filters

Available in Piston Base on request

