



Obturator Needle & Sling System

The Obturator Sling System is intended for the placement of a sub-urethral mesh for the treatment of female stress urinary incontinence (SUI) resulting from urethral hyper mobility and/or intrinsic sphincter deficiency



Product Details:

Monofilament Diameter	0.15 mm
Pore Size	1.3 x 1.0 mm
Burst Strength	162 psi

Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



Break Strength	MD 118 (ibs/in) CMD 111 (ibs/in)
Break elongation %	MD 107 CMD 175
Thickness	0.59 mm
Weight	100 gms/m ²

Description

Freedom Trans - Obturator Sling System

1 Set of Freedom Trans - Obturator Needles with CD

Kit of 1 Set Freedom Needles + 3 Freedom Trans - Obturator Sling System + CD



Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com