



Alak Oxygen Mask (Adult / Pediatric)



Features

- Made of clear, non-toxic, Food Grade PVC material
- Elastic strap provided to position the mask properly.
- Soft, large & designed in such a way that face mask fits comfortably.
- Specially designed for convenient oxygen therapy.
- Provided with 200 long star lumen Kink resistant tube to ensure continuous flow of oxygen.
- The tube has funnel shaped connector at the proximal end to facilitate simple & safe connection with the oxygen source.
- Sterile / Non Sterile packed in class 10000 clean room, Disposable / individual packed & Pyrogen free.
- It is designed to deliver oxygen to the patient during surgeries or critical care.