



Alak Ryle's Tube



Details

- The Tube is made of Excellent quality from non toxic PVC.
- Specially designed for naso gastric introduction for nutrition and aspiration of intestinal secretion.
- Closed distal end with metallic balls sealed into the tube to assist the passage of tube during intubations
- Four lateral eyes are provided for efficient aspiration and administration.
- Non irritant to delicate mucosa
- Tube with radio- opaque line for x ray visualization
- The tube is marked at 50 cm 60 cm &70 cm from the tip for accurate placement.
- Super smooth low friction surface facilitates easy intubation.
- Sterile / Disposable / individual packed.
- Note: For stomach tube length is 120 cms

Size: FG 8/10/12/14/16/18/20

Color Coded									
	Blue	Black	White	Green	Orange	Red	Yellow	Violet	Light Blue
Size FG	8	10	12	14	16	18	20	22	24