



## Skin Stapler

Skin Stapler which is simple designed and easy to operate. It is very reliable. No operational training is required. Staples with different external diameters and width can apply to sewing different operation incision. The sewing time is short and minimizes occurring of complication, and reduce operation cost. The range is sterilized with ETO.






### Product Details:

Type Wide	0.6mm Wire Diameter
Type Regular	0.55mm Wire Diameter
Span Regular	5.9 mm
Span Wide	7.2 mm
Leg Length Regular	3.9 mm
Leg Length Regular	4.9 mm
Packing - Prosec Skin Stapler	Box of 5
Packing - Staple Remover	Box of 20

### Application:

General Surgery, Heart and Bosom Surgery, Plastic Surgery, Burn Surgery, Neuron surgery, Bone Surgery, Gynecology and Obstetrics, Oncology Surgery, Field Rescue, Emergency Treatment



**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com