

SMART
Radiography X-Ray System
300mA/500mA





Wall Mount Panel

SMART RADIOGRAPHY X-RAY SYSTEM SPECIFICATION

MODEL NAME	SSL-2000	SSL-2000
POWER SUPPLY REQUIREMENT	1 Phase	1 Phase
INPUT RATING	400VAC/50/60Hz ± 10%	400VAC/50/60Hz ± 10%
LINE RESISTANCE	0.5 ohm	0.5 ohm
OUTPUT RATING	50 kVp	48 kVp
CONTROL PANEL	Microprocessor Controlled Pusher-Touch Operated	Microprocessor Controlled Pusher-Touch Operated
MAX. mA RANGE	40 - 150 mA	40 - 150 mA
MAX. mAs RANGE	2 - 200 mAs	2 - 200 mAs
MAX. mA RANGE	100, 1000, 2000, 3000, 4000	100, 1000, 2000, 3000, 4000
EXPOSURE TIME	0.02 - 60 Sec.	0.02 - 60 Sec.
X-Ray Tube	Rotatable Anode	Rotatable Anode
TUBEHOLD TYPE	Shel Type	Shel Type
RECTIFICATION	Full Wave Rectification	Full Wave Rectification
FOCAL SPOT	1.2 mm ² Small, 2.0 mm ² Large	1.2 mm ² Small, 2.0 mm ² Large
EXAMINEE FEET REST STAND	SSL-2F-00A	SSL-2F-00A
BUCKY TABLE	SSL-4 Way Floating Table	SSL-4 Way Floating Table
OPTIONAL*	SSL-4 Way Floating Table	SSL-4 Way Floating Table
	Strong Bucky 100-00d	Strong Bucky 100-00d
OPTIONAL BUCKY STAND	SSL-2F-00	SSL-2F-00
	Strong Bucky 100-00d	Strong Bucky 100-00d



Vertical
Bucky

Floor
Mount
Panel

- Optional -

- Table View -

