




Alak Stomach Tube



Features

- Specially designed to allow nasogastric introduction of nutrition and aspiration of intestinal secretion.
- Tube : Radio opaque, soft, kink free non toxic PVC tube
- Distal end is coned with four lateral eyes.
- Low friction Medical grade PVC reduces the incidence of nasal discharges.
- Proximal end is provided with universal funnel with shape connector.
- Marking: 50 cm, 60 cm, and 70 cm from the tip for accurate placement.
- Length : 125cm
- Sterile / Disposable / individual packed.

Sizes : FG 8 (Blue), 10 (Black), 12 (White), 14 (Green), 16 (Orange), 18 (Red), 20