



## Alak Suction Catheter with Thumb Control



### Details

- Suction Catheters® are suitable for removal of secretion from mouth, trachea and bronchial tubes.
- Manufactured from non-toxic non-irritant medical grade PVC
- Distal end is open with one lateral eye to avoid trauma during suction.
- Colour coded for instant easy size identification
- Frozen surface tubing for smooth intubation.
- Radio-opaque line,
- Provided with T type vacuum control valve with thumb control for proper maneuvering
- Available Plain (open/close tip), vacuum control valve with thumb control and syringe control facility.

Color Coded							
	Blue	Black	White	Green	Orange	Red	Yellow
Size FG	8	10	12	14	16	18	20