



Surgical Meshes for Hernia Repair

Polypropylene Meshes for Hernia Repair

Product Details:

COMPOSITION	POLYPROPYLENE		
CONSTRUCTION	MONOFILAMENT		
ABSORPTION TIME	NON-ABSORBABLE		
	<u>MESH 1</u>	<u>MESH 2</u>	<u>MESH 3</u>
MESH WIGHT	80G/M2	45G/M2	28G/M2
THICKNESS	0.50MM	0.43MM	0.30MM
PORE SIZE	0.5X0.7MM	1.3X1.0MM	1.0X1.0
ELONGATION	89%	77%	51%
MAXIMUM FORCE	25 KGF	19 KGF	15 KGF
SHELF LIFE	5 YEARS		

Registered Office: Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

Branch Office: 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

Tel: +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



ALAK

HEALTHCARE

AVAILABLE IN SIZE	6X11 CM, 7.6X15CM, 10X15 CM, 12X15 CM, 15X15 CM, 15X30 CM, 30X30 CM		
-------------------	---	--	--

