



Alak Twin Bore Nasal Oxygen Set Adult & Pediatric



Details

- Manufactured from Soft, non toxic PVC
- Designed for easy administration of oxygen and adjustable nasal tips for patient comfort.
- The mouth of the patient remains free for oral nutritional feeding and communication.
- Soft twin prong nasal tips to ensure equal volume of oxygen to both air passages.
- Star lumen tube ensures the supply of oxygen even if tube is accidentally kinked.
- Soft Funnel shaped connector facilitates easy connection to oxygen sources.
- Sterile/Disposable/ Individual Packed.