



**ALAK**  
HEALTHCARE

**ALDS-01 Plus**  
Digital Ultrasound System



- Multi-applications---Abdominal, Urology, GYN, OB, Pediatric, Small parts, Cardiology
- 10" CRT Monitor
- Dual Transducer connector
- THI (Tissue Harmonic Imaging)
- Weighs 10Kg
- Imaging Modes: B, B+B, B+M, M, 4B
- Two USB (2.0) ports, Video Port, VGA Port
- Automatic Diagnostic Report Systems

**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



## ALDS-03 Plus

### Digital B/W PWD Ultrasound System



- Multi-applications---Abdominal, Urology, GYN, OB, Pediatric, Small parts, Cardiology, Peripheral Vessels, Ortho
- 12.1" TFT LCD Monitor
- Dual Transducer connector
- THI (Tissue Harmonic Imaging) & Color Chroma
- Optional Battery back up
- Weighs 6 Kgs
- Pulse wave Doppler available with Doppler calculations
- Imaging Modes: B, B+B, B+M, M, 4B, PWD
- Image Storage in Local Disk & USB
- 256Frames Cine loop



## ALDS-05 Plus Digital Hand Carried Color Doppler



- Imaging mode: B, B+B, 4B, M, Color, PWD, HPRF, Power, Dir Power Doppler, DUAL LIVE (B & Color) & Duplex / Triplex Imaging
- 12.1 inch LCD Monitor
- 2 Transducer connector
- Speckle Reduction Algorithm(SRA)
- THI, TSI, digital beam forming
- Spatial Compound Imaging
- Phase inversion Harmonic Imaging
- Multi frequency Transducers
- Software packages: OB, GYN, Cardiac, Small Parts, URO, Peripheral, Carotid, Thyroid, MSK, Vascular
- Standard battery back up
- 256 frames CINE loop memory, 2 USB ports

**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com





## ALDS-05 Elite Portable Color Doppler



- 15" LED Monitor
- Advanced Portable Color Doppler
- Backlight Control Panel
- Standard Dual Transducer Connector
- Various Storage & Printing Possibilities
- Standard 2P Connector
- Internal 500GB HDD
- Built-in Battery (2 Hrs Backup)
- Measurement and Calc. Package
- Direct Multiple Image Print out



## ALDS-015 Plus Digital Color Doppler



- Imaging mode: B, B/B, B/M, 4B, Color Doppler Flow Imaging Pulse Wave Doppler, Power Doppler Imaging, Directional Power Doppler
- 15" inch LCD Monitor
- 4 Ports (2 Active ports & 2 Parking ports)
- Speckle Resistance
- Multi-beam Forming, Beam forming optimization
- Phased Inversion HCI and Spatial Compound Imaging
- Large Capacity Integrated Hard Disk
- Multi frequency Transducers
- Compound Technology
- Software packages: ABD, OB, GYN, Cardiac, Small Parts, URO, Vascular, Anesthesia

**Registered Office:** Survey No. 81, Plot No. 16-17, Junagadh Road, Near Lion's School, Jetpur - 360370, Gujarat, India.

**Branch Office:** 401, Vithal Bhavan, 1st Floor, Chira Bazaar J.S.S. Road, Mumbai – 400 002, India.

**Tel:** +91 022 66370811, +91 022 66370812 **Email:** info@alakhealthcare.com **Website:** www.alakhealthcare.com



**ALAK**  
HEALTHCARE

**Ultrasound Scanner for Veterinary also available**

**DS01Plus Vet**



**DS03Plus Vet**



**DS05Plus Vet**

